White Novelys [®] SNEAKER CARE RANGE

TECHNICAL CHARACTERISTICS

Repigments, renovates and beautifies.

- Fluid, easy to apply.
- Highly covering, thanks to a very high concentration of pigments.
- Long-lasting resistance to rubbing and water contact.

CARE INSTRUCTIONS

1. Brush the article to remove any dust or dirt

- 2. Shake the bottle to homogenize the product
- 3. Apply the sponge upside down, pressing lightly
- 4. Dry for 15 minutes

5. Rinse the applicator sponge under water and cap the bottle.

If necessary, clean smooth white leather with RENOMAT (Ref.0514) Make a test on a hidden part beforehand. Rinse the foam from the applicator to preserve it. Waterproof treated articles with INVULNER aerosol (Ref.0745)

For other smooth white leather items, to be recolored or renovated, use : - WHITE RENOVATOR CREAM (Ref.085121)

- white JUVACUIR (Ref.080321)

WRAPPING / PACKAGING

VOLUME	WRAPPING			PACKAGING			
	TYPE	DIMENSION IN MM		QUANTITY PER	DIMENSION IN MM		
		DIAMETRE	HIGHT	PACKAGE	LENGHT	WIDHT	HIGHT
75 ML	APPLICATOR BOTTLE	35	135	12	110	75	150

NEAKE

CARE

ARTICLE CODE / GENCOD

DESIGNATION	VOLUME	GENCOD			
DESIGNATION		ARTICLE CODE		E CODE	
				TARIF CODE	
WHITE NOVELYS	WHITE	75 ML	332401	0303	218







Tarif code: 0303

WHITE

NOVELYS

CUIR, TEXTILE BLANC LEATHER, SYNTHETIC WEISS, WIT, BLANCO

NET WT. 3.7 oz (105g)

New

VOLUME :

75 ml

Data Sheet